

For highly sensitive detection

# Gold Nanoplates


## Features

- 1) An aqueous dispersion of disk shaped Gold Nanoparticles.
- 2) Exhibits a vivid blue color due to absorption caused by localized surface plasmon resonance.
- 3) By applying Gold Nanoplates to diagnostic agents, it is expected to enable highly sensitive detection and reduce the cost of antibodies.


## Characteristics of Gold Nanoplates

### Shape

Gold Nanoplates are disk shaped triangle nano particle



Electron microscope image




### Absorption characteristics

Gold Nanoplates have a sharp optical absorption waveform and exhibits a vivid blue color. The maximum absorption wavelength of Gold Nanoplates are determined by the particle size and aspect ratio (maximum length/thickness).



Absorption spectrum



Gold Nanoplates	
Name	Au-WPPLC1-C
Dispersion medium	Water
Dispersion materials	Trisodium citrate, etc.
Gold content	0.04~0.06 mg/g
Absorption properties (The maximum absorption wavelength)	615 ± 15nm
Particle size (Maximum length)	Approximately 50 nm
pH <sup>※1</sup>	4~7

※1 pH: HORIBA, Ltd. Twin pH meter

※2 Electron microscope image: Hitachi, Ltd. Scanning Electron Microscope SU 70

※3 Spectral characteristics : Diluted with water to an arbitrary concentration and measured. (Shimadzu Corporation UV Visible Spectrophotometer MPC3100UV)

Note: The above values are for reference only and are not standard values.


For highly sensitive detection

# Gold Nanoplates

## Example of the specimen detection

### 1, Immunochromatography method

By combining Gold Nanoplates with antibodies, The complex can be used as a coloring agent in immunochromatographic tests. The detection line shows a clear color tone of the dispersion liquid. Also, the same sensitivity can be obtained with a smaller amount of antibody compared to spherical gold nanoparticles.




Human chorionic gonadotropin (hCG) detection result (Visual judgement)

Metal nanoparticles	Antibody addition amount [μg/mL(Abs. 1.0)]	Amount of hCG [mIU]			
		100	10	1	Blank
① Gold Nanoplates (Au-WPPLC1-C)	1.25	++	++	+	—
② Spherical gold nanoparticles (Commercial product)	2.50	++	++	+	—


Judgement results : ++...Strong positive + ...Positive --- Negative

hCG detection result (immunochromatographic test strips after test)




### 2, Agglutination method

The antigen can be detected by observing the change in absorption characteristics of nanoplate-antibody complex, because of the complex interact with the antigen. Gold nanoplates show a more conspicuous change than spherical gold nanoparticles.



Schematic diagram of the agglutination method

The nanoplate-antibody complex agglutinate around the antigen and precipitates.



Changes in spectral characteristics over time after addition of hCG (Inset: color change)

When the nanoplate-antibody complex agglutinates and precipitates due to the antigen antibody reaction, two changes occur.

- ① Decreased extinction of maximum absorption wavelength.
- ② Increased extinction in the long wavelength range.